



FotonaSmooth™

First Non- surgical Laser for Gynaecologists

- **IntimaLase™** - laser vaginal tightening
- **IncontiLase™** - stress urinary incontinence laser treatment
- Minimally invasive, incisionless
- No anaesthesia required
- Walk-in/walk-out procedure
- **FotonaSmooth™** technology that will revolutionize gynaecology forever

Fotona[•]
choose perfection

The Highest Performance, Best Made Laser Systems in the World

A true incisionless laser treatment for vaginal relaxation syndrome

IntimaLase™

- laser vaginal tightening (LVT) through photothermal tightening of the vaginal canal
- minimally invasive, incisionless
- no anesthesia needed
- safe, quick and easy ambulatory procedure
- walk-in/walk-out procedure

How does IntimaLase™ work?

- Fotona's 2940 nm Er:YAG non-ablative laser with proprietary "SMOOTH-mode" technology
- Thermal effect on the vaginal mucosa that stimulates collagen remodeling and the synthesis of new collagen fibers in the vaginal mucosa tissue and collagen-rich endopelvic fascia.
- The final result of collagen neogenesis and remodeling is tightening of the vaginal canal.



Vaginal relaxation syndrome



Tight vaginal canal

A minimally invasive solution for mild and moderate stress and mixed urinary incontinence

IncontiLase™

- laser treatment of mild and moderate stress and mixed urinary incontinence (SUI)
- photothermal tightening of the urethral and anterior vaginal wall region
- minimally invasive, incisionless
- no anesthesia needed
- safe, quick and easy ambulatory procedure
- walk-in/walk-out procedure

How does IncontiLase™ work?

- Fotona's 2940 nm Er:YAG non-ablative laser with proprietary "SMOOTH-mode" technology
- Collagen remodeling and the synthesis of new collagen fibers in the region of anterior vaginal wall as well as the vestibule and urethral orifice due to the thermal effect on the vaginal mucosa.
- Greater support to the bladder and the return of normal continence function through collagen neogenesis and shrinking and tightening of vaginal mucosa tissue and collagen-rich endopelvic fascia.



Fotona SMOOTH mode treatment of the anterior vaginal wall



Mild and moderate stress and mixed urinary incontinence



After IncontiLase treatment